



## Hybridization Solution II Product Specification Sheet

**CAT. NO.:** IH-60010 (HYB-0002)

**Description:**

Hybridization Solution II may be used for pre-hybridizations and hybridizations of nucleic acid probes to nucleic acids immobilized on nitrocellulose or nylon membrane. UP to 50% deionized formamide may be added for RNA hybridizations

**Quantity:** 200ml.

**Formulation:**

20% Formamide, 3X SSC, 5X Denhardt's.

**Storage:**

It is stable at least one year when store at 4-8°C.

**Note:**

The product is for research use only, not for human or drug use.

# Maxim Biotech, Inc.

780 Dubuque Avenue, South San Francisco, CA 94080 USA  
Tel: 650-871-1919 Fax: 650-871-2857 E-mail: [mbi@maximbio.com](mailto:mbi@maximbio.com)